

STEPS
Sunrise Training & Education Programs

Dependent Mobility Bases: How do I choose?

Presented by:
Angie Kiger, M.Ed., CTRS, ATP/SMS

MANUAL
POWER
ADULT
PEDIATRICS
SEATING
GERIATRICS
CONTROLS
FUNDING

Objectives

Upon completion of this presentation participants will be able to:

- List 3 indicators as to when it is appropriate to recommend an adaptive stroller.
- List 3 clinical indicators for use of a tilt-in-space wheelchair instead of an adaptive stroller.
- List 3 seating options available on a tilt-in-space wheelchair.

STEPS
Sunrise Training & Education Programs

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

2

Sooo....What's the Big Deal?

<p>Proper Seating Impacts:</p> <ul style="list-style-type: none"> - PLAY! - Access to environment - Developmental milestones - Reduction or prevention of risk for injuries in the future - Respiratory status - Gastrointestinal function 	<p>Mobility Impacts:</p> <ul style="list-style-type: none"> - PLAY! - Cognition - Decision making - Social interaction - Spatial awareness - Body control in gravity
---	---

STEPS
Sunrise Training & Education Programs

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

3

What's The Big Deal?

“The first time I saw her sitting in her KidKart I was shocked at how upright she was sitting and how much more engaged she was in her surroundings. It was a wake-up call to her dad and I about how important proper positioning is for Maryn.”

- Kim (Arlington, VA)

“Providing clients with proper seating and positioning is vital for feeding, communication, and socialization!”

- Becky, CCC-SLP (Pensacola, FL)

STEPS
Sunrise Training & Education Programs

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

4

Perspective

“I understand why my parents want to see me walk and talk, but it's exhausting to live in a body that feels like the property of everyone else.”

- Martin Pistorius

STEPS
Sunrise Training & Education Programs

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

5

Slicing Up The “Pie”

Who is fighting for a piece of the pie when it comes to choosing equipment?

STEPS
Sunrise Training & Education Programs

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

6

Setting Goals

❖ Targeted areas to consider:

- Positioning
- Stability
- Manage Tone
- Changes in size
- Sitting tolerance
- Skin integrity
- Function

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Special Training & Educational Programs

7

Goals to Keep In Mind

OTHER AREAS OF CONSIDERATION

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Special Training & Educational Programs

8

What are the goals?

- **Child**
 - Peer interaction
 - Independence
 - Fun
 - Play
 - Explore
 - Interact
 - Learn - feel, touch, do
 - Looks "cool"
- **Family**
 - Aesthetics (low profile)
 - Acceptability
 - Accessibility
 - Ease of use
 - Comfort
- **Clinician**
 - Good positioning
 - Complimenting therapy goals
 - Easy to use
 - Promote independence
 - Safety
- **Funding Source**
 - Thorough documentation
 - Meeting the criteria




CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Special Training & Educational Programs

9

Question...



Is there EVER conflict???

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Special Training & Educational Programs

10

Recommended Flow

- Explain the purpose
- Review goals
- Complete history
 - Medical
 - Dx
 - Rx
 - Surgeries
 - School
 - Therapy
 - Equipment & AT experience
 - Environment & transportation
 - Likes/Dislikes
- Observation
 - Before introductions
 - Current equipment
 - Out of equipment
- Hands-on
 - Skin
 - Mat evaluation
 - Cognitive skills
 - Sensory status
 - Function
- Featuring matching
- Education
 - Funding process
 - Resources for support
 - Planning for lifetime of equipment
- Recommendations
 - Appointment summary
 - Contact information
 - Follow-up plan
 - Review goals...again!

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Special Training & Educational Programs

11

Additional Evaluation Considerations

- **Screenings, reports, or formal consultation**
 - Physicians – neurologist, ophthalmologist, orthopedic surgeon, physiatrist, etc.
 - Speech Therapist
 - Teacher - classroom aide
 - School therapists – IEP
 - Audiologist
- **Cognition**
 - Ability to following directions (one-step, two-step, multi-step, related, non-related, etc.)
 - Initiation of exploring the environment independently
 - Visual learner vs. auditory learner (supports needed?)
- **Auditory status**
 - Acuity vs. processing
- **Communication Status**
 - Verbal vs. non-verbal
 - Picture based vs. word based

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Special Training & Educational Programs

12

Additional Evaluation Considerations

- **Vision Status**
 - Acuity vs. processing
 - Field loss, field neglect, color blind, visual motor, etc.
- **Cortical Visual Impairment**
 - CVI is a neurological condition that is the leading cause of visual impairment of children in the US.
 - Commonly behaviors associated with CVI include:
 - Require movement to see
 - Visual field preference—see better looking at objects in certain directions (such as left or right periphery).
 - Lack of visual-motor match—look and touch occur as separate functions (looks, turns head away from item, then reaches for it)
 - Light gazing and/or non purposeful gaze—often need more light to see or will gaze non-purposefully because they can't make sense of what they see

CVI Common Comorbidities

- Cerebral Palsy
- Cognitive impairment
- Seizure disorder
- Microcephaly
- Hearing loss
- Memory dysfunction
- Hyper or hyposensitivity to sensory stimulation



Time For The Stroller Talk

- When is it appropriate to recommend adaptive seating for a child?
- What is the best way to approach the family?
- What information do you need to have prepared prior to talking with the family?
- Why might a family tell you “no”?

Standard Strollers



Early Intervention Adaptive Strollers



Standard vs. Adaptive

- | <u>Standard Stroller</u> | <u>Adaptive Stroller</u> |
|--|---|
| <ul style="list-style-type: none"> • Mainstream • Lower profile • Standard seating option only • Limited recline available in some brands. • Easily transportable | <ul style="list-style-type: none"> • Seating & positioning options. • Durability • Storage and accessories for transporting medical equipment • Seating system removable • Tilt, recline, & transit available • Foldable base |



Seating Considerations

- Seating
 - Reversible?
 - Size?
 - Growth?
 - Weight limits?
 - Cushion, back, positioning components, headrest options?
 - Tilt, recline, ELRs?



Frame Considerations

- Transit
- Folding
- Weight
- Age appropriate



Accessories



Stroller To Wheelchair

- Why do some parents hesitate with transitioning from a stroller to a wheelchair?
 - Strollers (even adaptive strollers) look more mainstream
 - Accessibility
 - Home
 - Transportation
 - Ease of getting from point A to point B
 - Acceptance
 - Funding
- Why is it important to move a child from a stroller to a wheelchair?
 - Positioning
 - Age appropriateness
 - Seating
 - Access to the environment

Tilt-In-Space Wheelchairs



Points to Consider

- Seating compatibility
 - Standard
 - Custom
 - Tilt
 - Recline
- Accessories
 - Headrests
 - Upper extremity supports
 - Lower extremity supports
 - Positioning belts
 - Wheel options
- Weight
- Folding options
- Configuration
- Growth



Tilt Options

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Research Program

25

Folding Options

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Research Program

26

Consider A Seating Ladder

Custom made seating – molded

Custom made seating – linear and contoured

Off the shelf seating – customizable

Basic off the shelf seating – non customizable

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Research Program

27

Off the Shelf Cushions - Customization

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Research Program

28

Cover Considerations

- Material
- Stretch
- Layers
- Style

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Research Program

29

Maintenance & Care

- Concerns about moisture or incontinence?
 - Type of protection needed
- How is the cushion washed/cleaned?
 - Additional covers?

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Research Program

30

Out of the Box Backs & Adaptations

STEPS
Seating Training & Assessment Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

31

Custom Fabricated Seating: Unique Postures

STEPS
Seating Training & Assessment Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

32

Custom Seating Surfaces

STEPS
Seating Training & Assessment Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

33

Custom Backrest Shapes

STEPS
Seating Training & Assessment Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

34

Foam In Place

- Back rest/Carrying Structure
- 1" Sunmate foam
- Loose 4" deep cover
- 1 Unit of liquid soft Sunmate foam.
- Common Mod's
 - Lateral Supports
 - I style back
 - Curved wood
 - Custom base foam

STEPS
Seating Training & Assessment Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

35

"Hybrid" Option - ShurShape

- What is it? What is the process?
 - Chopped up foam pieces in a large vacuum sealed bag (one for the seat and one for the back).
 - Glue is poured in the bag to mix with the foam. Next a suction machine is hooked up to the bag and the glue/foam concoction is molded to the client.
 - A "dry run" can be done by using the suction without the glue to determine how the system will impact the client.
- Things to keep in mind
 - Slight adjustments can be made while mold is setting.
 - Carrying structures enable specific seat depth, back height and angle adjustments to accommodate an optimum fit to the wheelchair base.
 - Can do a dry mold
 - The foam can be cut away or modified in the field by a seating professional.
 - Mixture must be mixed thoroughly.
 - Not ideal for clients with high tone.

STEPS
Seating Training & Assessment Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

36

Custom Molded Seating

- In **some** cases it is easier to attain consistent positioning/targeting with a mold needing less adjustment or repositioning after transfer.
- Molding systems allow direct involvement of the clinician in the final shape as opposed to taking measurements.
- Molds eliminate moveable, adjustable parts – **can be positive**
- Molds eliminate adjustment and fine tuning – **can be negative**



CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

37

Additional Seating Components



CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

38

Lateral Thoracic & Pelvic Pads



Thoracic



Curved



Pelvic



Tapered



CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

39

Lateral Mounting Hardware



Fixed



Adjustable



Swing-away



Removable




Hardware cover



CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

40


Upper Extremity Supports



CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

41

Medial Thigh Supports and Hardware



Wedge



Oval



CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

42

Lower Extremity Supports

Contracture Footrest System
Inside Mount Outside Mount

STEPS
Special Training & Research Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

43

Lower Extremity Positioning Supports

Shoe holders Leg Extension

Foot Boxes Knee Adductor Stabilizer

STEPS
Special Training & Research Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

44

Positioning Belts

STEPS
Special Training & Research Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

45

Headrest Pad & Component Options

STEPS
Special Training & Research Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

46

Dynamic Seating Components

Whitmyer Flex Interface Bracket

Millers Dynamic Footrests

STEPS
Special Training & Research Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

47

Growth Considerations for Mobility Bases

- What is too much growth?
- Chair frames
 - Should have both width, depth and seat to floor height change capability
- Quick Adjustments
 - Moving the back post
 - Growing the cross brace, or strut tubes
 - Swapping side frames (height)
- More Involved Adjustments
 - Adjustment kit
 - New frame

STEPS
Special Training & Research Program

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

48

WAIT... Don't Forget The Most Important Parts!

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Assessment Program

49

Dependent Transport Strollers

- Frame Considerations
 - Transit
 - Folding
 - Weight
 - Length of time

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Assessment Program

50

Dependent Transport Strollers

- Seating Considerations
 - Positioning support
 - Adjustments
 - Tilt
 - Recline

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Assessment Program

51

Critical Questions

- Who is the funding source?
- What is the client's medical history?
 - Diagnosis (primary, secondary, etc.)
 - Surgeries (previous and upcoming)
 - Medications (past, present, future)
- What equipment has the patient had?
 - Not just wheelchairs
 - When was it received, why does it no longer meet their needs (medical - primary)? Who funded the equipment?

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Assessment Program

52

Documentation

- Is your evaluation completed electronically or hand-written?
- Letter of medical necessity?
- Template or not?
- Important documentation reminders:
 - Your clients are individuals
 - Proof-read!
 - Contradictions

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Assessment Program

53

Final Questions to Ask Yourself

1. Does the client and/or caregiver have a clear understanding of the plan?
2. Have I specified that the recommended equipment is in fact the minimal equipment essential to this client?
3. Have I demonstrated how I ruled out lesser level equipment?
4. Is the equipment that I am recommending in fact the least costly alternative?
5. Do I have all of the information needed for funding?
6. Has my documentation left the reader with a clear picture of the consequences to the client in the absence of having the recommended equipment?

CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/15

STEPS
Senior Training & Assessment Program

54

Words of Wisdom

"Can't is not a bad word. We ALL have "can't's" but when there is a can't that significantly impairs your life, you find a **work around**. So for the little boy who couldn't walk a wheelchair is simply a **work around**."

Dr. Lisa Thornton,
Pediatric and Adolescent Physiatrist
LaRabida Children's Hospital & Schwab Rehabilitation Hospital



CONFIDENTIAL AND PROPRIETARY Duplication or Distribution Prohibited 11/12/13

55



Contact Information

[@ATigerKiger](https://twitter.com/ATigerKiger)
angie.kiger@sunmed.com

www.EducationInMotionBlog.com

For more information please visit our website at: www.numotion.com

Course transcripts and additional copies of certificates of completion are available upon written request.
5501 Wilshire Blvd NE, Ste C, Albuquerque, NM 87113 1.800.500.9150.

