

# THE CASE FOR WHEELS: Selecting the Clinically Appropriate Power Wheelchair Base

Julie Piriano, PT, ATP/SMS  
Quantum Rehab  
800-800-8586  
[jpiriano@quantumrehab.com](mailto:jpiriano@quantumrehab.com)

---

---

---

---

---

---

---

---

## Agenda

- Perceptions of Drive Wheel Position
- Advantages/Disadvantages of FWD, MWD and RWD
- Suspension
- Tracking Technology
- Making an Informed, Individual Choice

---

---

---

---

---

---

---

---

## Learning Objectives

Participants will be able to:

- List two advantages of each of the 3 main power wheelchair configurations available.
- Differentiate the performance of front, mid and rear wheel drive wheel configuration across "outdoor terrain".
- Relate a minimum of two reasons why suspension would be medically necessary for an individual with a neurological impairment.
- Justify the medical need for tracking technology on a power wheelchair with a specialty control device.

---

---

---

---

---

---

---

---

### Satisfactory Completion Requirements

- Individual course registration via GoSignMeUp (GSMU) account
- 100% in-class attendance
- Once instructor finalizes attendance, you will receive an email to complete a post-course survey
- Upon completion of survey, your certificate is immediately emailed to you or available in your GSMU account
- No partial credit will be awarded



---

---

---

---

---


---

---

---

### Audience Participation

- WiFi –
- Go to [www.menti.com](http://www.menti.com)
- Once at that website there will be a box to fill in a code to participate
  - The code required is **44 07 79**
  - This will bring you to the questions for today
  - You will be allotted 5 seconds to answer each question by selecting your answer and hitting **SUBMIT**



---

---

---

---

---

---

---

---

### Question 1

What drive wheel configuration is the best overall for individuals with disabilities?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



© mark da tait. www.markda.com

---

---

---

---

---

---

---

---



**Pride**  
Learning Institute

Is There a Drive Wheel Configuration that is "Best" for an Individual's Life Style?

---

---

---


---

---

---

---

---



**Pride**  
Learning Institute

There is no one mobility base or seating system configuration that is appropriate for an individual with a mobility impairment based on diagnosis, prognosis, deformity, disability or level of function.

**It is as unique as the individual!**

---

---

---

---

---

---

---

---



**Pride**  
Learning Institute

**Life is NOT a Level Tile Floor!**

- When, where and how will they will use the power wheelchair?
  - Need to consider **ALL** settings of anticipated use!
- What features are going to be most important?

---

---

---

---

---

---

---

---

## In the Beginning...

- There was only rear wheel drive...
- Then front wheel drive was introduced
  - Initially it didn't handle as well as RWD in places where active RWD users were use to taking their chairs but...
  - It handled indoor mobility, turning and maneuverability better than RWD = **big swing toward FWD recommendation**
- Mid Wheel drive was introduced last
  - Resolved many issues found by power wheelchair users in a wide array of settings
  - It was stable AND maneuverable = **shift towards MWD recommendation**
- As a result the RWD configuration was almost extinct
- Many FWD prescribers were hesitant to shift to MWD recommendation due to familiarity with FWD configuration



---

---

---

---

---

---

---

---

## Today's Power Wheelchairs!

- What are the deciding factors?
- How do we inform and educate our consumers, and **empower** them to make the "correct" or most advantageous choice?
- All they see is a **BLUE** chair!
- Are WE informed, educated and **OBSERVANT** of their needs and capacities?



---

---

---

---

---

---


---

---

## Question 2

What drive wheel configuration is most maneuverable in all environments?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---


---

---

### Question 3

What drive wheel configuration is best in tight spaces?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---

---

---

### Question 4

What drive wheel configuration is best out doors?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---

---

---

### Question 5

What drive wheel configuration is best for tracking straight?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---

---

---

### Question 6

What drive wheel configuration is best at high speeds?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---

---

---

### Question 7

What drive wheel configuration handles rougher terrain best?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---


---

---

### Question 8

What drive wheel configuration is most responsive when driving forward into a turn?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---

---

---

### Question 9

What drive wheel configuration performs best going uphill?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---

---

---

### Question 10

Which drive wheel configuration climbs obstacles best?

- A. Rear wheel drive
- B. Mid/center wheel drive
- C. Front wheel drive



---

---

---

---

---

---

---

---

What do you think are the advantages of front wheel drive?



---

---

---

---

---

---

---

---

### Advantages of Front Wheel Drive



- Make tight turns around corners
- Climbs obstacles well
- Large front wheels drive over obstacles with less impact on the user
- As long as the front anti-tip wheels are high enough FWD will transition over grade changes easily
- Provides an advantage for individuals with tight hamstrings for lower extremity positioning
- Works well for people who's center of gravity tends to be "forward" in the seating system. (i.e., bariatric population)
- Because FWD pulls it handles softer terrains better (i.e., grass, gravel, etc.)

---

---

---

---

---

---

---

---



What do you think are the disadvantages of front wheel drive?

---

---

---

---

---

---

---

---

### Disadvantages of Front Wheel Drive

- May require slower driving speeds to maintain control
- Anti-tip casters can interfere with foot placement for stand pivot transfers
- Turning in tight spaces can be tricky for new users because more chair is behind them
- Turning quickly can cause dizziness because the inner ear is further from the drive axle
- Harder to drive straight at higher speeds
- Has a tendency to veer uphill traversing side slopes
- Less control with non-proportional controls




---

---

---

---

---

---

---

---



What do you think are the advantages of mid wheel drive?



---

---

---

---

---

---

---

---

### Advantages of Mid-Wheel Drive



- Tightest turning radius for a 360° turn
- Most intuitive to drive for individuals who have previously been ambulatory
- Offers comparable control to rear wheel drive
- Climbs obstacles fairly well, can be limited by how high castor wheels lift
- Having 6 wheels on the ground provides stability to the base
- As long as the front caster suspension has enough travel the MWD w/c will transition over grade changes easily
- Good traction on most surfaces, inclines and side slopes

---

---

---

---

---

---

---

---



What do you think are the disadvantages of mid wheel drive?

---

---

---

---

---

---

---

---

### Disadvantages of Mid-Wheel Drive

- Front caster wheels can interfere with stand pivot transfers
- Front caster wheels can be problematic for individuals with tight hamstrings
- Can high center on very uneven terrains
- Have limitations in the height of obstacle they can climb
- Can give a sense of "pitching forward"
- Jolting forces directly blow the user



---

---

---

---

---

---

---

---

What do you think are the advantages of rear wheel drive?



---

---

---

---

---

---

---

---

### Advantages of Rear Wheel Drive

- Provides very good control at higher speeds
- Handles aggressive terrain very well
- Largest base = Most stable
- Will transition grade changes easily
- Can allow a care provider to assist with climbing a curb or obstacle due to drive wheel location in rear
- Can provide enhanced control with non-proportional inputs
- May be familiar for some users



---

---

---

---

---

---

---

---



What do you think are the disadvantages of rear wheel drive?

---

---

---

---

---

---

---

---

**Disadvantages of Rear Wheel Drive**

- Users typically have greater access issues due to large turning radius
- Front caster wheels can interfere with footplate placement
- May require some speed to get over obstacles that must be climbed
- Must drive straight on, to climb higher obstacles
- Front caster wheels can be problematic for individuals with tight hamstrings
- Least intuitive to drive for new users who were previously ambulatory
- Tend to veer downward traversing a side slope

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

Ⓞ Geometric Centre  
 Ⓞ Centre of Gravity  
 Ⓞ Centre of Head  
 Ⓞ Centre of Rotation

Pride Learning Institute

---

---

---

---

---

---

---

---

### Turning Circles

RWD      CWD      FWD

Pride Learning Institute

---

---

---

---

---

---

---

---

### Shape of Space

RWD      CWD      FWD  
 Low ratio    High ratio    High ratio    Low ratio

- Generally a CWD or HR chair provides a reasonable compromise with tight space requirements for effective maneuverability

Wheels counter rotating

One wheel fixed

Pride Learning Institute

---

---

---


---

---

---

---

---




---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

### Which is Best Suited for Indoors vs. Outdoors?

	Indoor		Outdoor	
	AB	ID	AB	ID
LR RWD	Red	Yellow	Red	Yellow
HR RWD	Red	Green	Red	Green
CWD/MWD	Green	Green	Yellow	Green
HR FWD	Yellow	Green	Green	Yellow
LR FWD	Yellow	Yellow	Green	Red

■ Least advantageous   
 ■ Acceptable performance   
 ■ Best performance

AB = Allan Boyd, (2011) Rehab Magazine Comparing PWC Platforms  
 ID = Ian Dennison, (2012) Understanding the Potential of Powered Mobility Devices, Proceedings of the 2012 International Seating Symposium, PPT presentation

---

---

---

---

---

---

---

---

### Performance Analysis – Indoor

- Speed
- Handling
- Size
- Maneuverability
  - Rate of turn
  - Turning radius
  - 3 point turn
  - Shape of space
  - Corridors and doors
  - Visual cues and ease of learning

	Indoor	
	AB	ID
LR RWD	Least advantageous	Acceptable performance
HR RWD	Least advantageous	Best performance
CWD/MWD	Best performance	Best performance
HR FWD	Acceptable performance	Best performance
LR FWD	Acceptable performance	Acceptable performance

■ Least advantageous 
 ■ Acceptable performance 
 ■ Best performance

---

---

---

---

---

---

---

---

---

---

### Performance Analysis - Outdoor

- Speed
- Handling
- Curbs / Curb cuts
- Acceleration / Braking
- Tracking
- Changing Direction
- Range
- Soft Terrain
- Comfort

	Outdoor	
	AB	ID
LR RWD	Least advantageous	Acceptable performance
HR RWD	Least advantageous	Best performance
CWD/MWD	Acceptable performance	Best performance
HR FWD	Best performance	Acceptable performance
LR FWD	Best performance	Least advantageous

■ Least advantageous 
 ■ Acceptable performance 
 ■ Best performance

---

---

---

---

---

---

---

---

---

---

### Question 11

Does suspension make a difference when considering drive wheel configuration?

- A. Yes
- B. No




---

---

---

---

---

---

---

---

---

---

### What does Suspension Provide?

- Absorption of jolting /vibratory forces on wheelchair and wheelchair user
  - Influences an individual's posture and balance in the (stability of positioning)
  - Reduces or mitigates spasticity and/or reflex activity
  - Supports pain management
  - Increases sitting tolerance
  - Improves durability of the power base
- Ability to maneuver over a variety of terrains, optimizing environmental transitions
- Stability of the power base



---

---

---

---

---

---

---

---

### Question 12

What industry did Quantum consult with to help design the suspension system in the 4FRONT and Edge 3?

- A. Mountain bike
- B. Motocross
- C. Construction vehicle
- D. Monster truck



---

---

---

---

---

---

---

---



---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

### Impact of Suspension on the Individual

- “Vibration, shock (single event and repeated), and motion have a significant effect on the health and quality of life for individuals who utilize a wheelchair for mobility.
- Vibration and shock can cause back pain and injury, which has been well documented in the literature.”



(Griffin, 1975, 1990; Kitazaki & Griffin, 1998; Paschold & Mayton, 2011; Wilder, Magnusson, Fenwick, & Pope, 1994)

---

---

---

---

---

---

---

---

### Impact of Suspension on the Individual

TRUE shock absorption limits the jolting/vibratory forces, which decreases the terrains effect on the individual:

- Support postural control, maintain positioning in the chair and contact with support components
- Preserve access to the drive control device, especially a specialty control
- Manage spasticity
- Minimize pain
- Decrease fatigue and increase sitting tolerance
- Reduce fear of driving across rough terrain and increase access to all environments of use



---

---

---

---


---

---

---

---





### Position of the Person in the Wheelchair – Center of Gravity

- Center of gravity of the individual in the seating system impacts the performance of the chair
- The center of gravity will significantly alter the performance of the suspension of the wheelchair
- If casters in the front or rear of a wheelchair are overloaded it affects wheelchair performance and drivability
- Front or rear loading of the caster wheels will affect the durability of the tires/wheels themselves.
- Note: People who carry a great deal of weight forward may get the best performance out of a front wheel drive

---

---

---



---

---

---

---

---

### Tracking Technology – What is it?

A **feature** on the motors that interface with the control module so that the speed of the drive wheels can change in relation to one another to:

- Maintain a straight course of operation
- Compensate for discrepancies in terrain, transitions and thresholds
- Track independently over surfaces the consumer encounters in their home and community

---

---

---



---

---

---

---

---

### Tracking Technology – What does it do?

- Allows an individual to control their w/c without having to make continual adjustments due to unequal drive wheel revolutions
- Minimized the need for excessive movements of the (hand, head, neck, chin, foot or other body part) to stay on track
- Compensates for veering when using non-proportional (switch) type controls
- Promotes safe, accurate operation of the w/c

---

---

---

---

---

---

---

---



**Pride**  
Learning Institute

### Tracking Technology – Why is it Medically Necessary?

- Allows individuals with limited motor control, muscular fatigue, physiological fatigue, balance challenges and/or spasticity to compensate for their loss of function
- Increases length of time individuals can drive
- Minimizes the risk for accidents or an adverse occurrence

12

---

---

---

---

---

---

---

---



**Pride**  
Learning Institute

After EDUCATION and TRIAL your Client Needs to Decide which is BEST for THEM!

13

---

---

---

---

---

---

---

---

### ASK Questions, TEST Drive Chairs, Know WHAT you are recommending and WHY!

**Pride**  
Learning Institute

**The Case for Wheels Questions**

1. What drive wheel configurations does each manufacturer offer? (circle all that apply)

AVID	Invacare	Permobil	Quantum	Rovi	Sunrise
FWD /MWD /RWD	FWD /MWD /RWD	FWD /MWD /RWD	FWD /MWD /RWD	FWD /MWD /RWD	FWD /MWD /RWD

Notes

14

---

---

---

---

---

---

---

---

## Thinking Out Loud

- What are two things you learned and how can you relate it to your practice?
- How does this course relate to your practice setting?
- How will this course change your behavior in the future?

---

---

---

---

---

---

---



QUANTUM REHAB®

## Thank You for Your Time – Any Questions?

Quantum Rehab  
Clinical Education Department  
E: [education@pridemobility.com](mailto:education@pridemobility.com)

---

---

---

---

---

---

---