

Leading the way toward independence in pediatric mobility

Continuing Education Completion Requirements for 0.10 CEU (1-hour)

Attendees must: Attend the entire course • Sign In/Out • Complete the assessment within 2 weeks Assessment(permobil.com)

Permobil Representative Locator https://hub.permobil.com/sales-locator


1) Motor Milestones¹ (birth- 2 years old)

*CDC's 2022 guidelines do not include a crawling milestone.









2 months Holds head up when on tummy, moves both arms and both legs, opens hands briefly

4 months Holds head steady without support when you are holding her, holds a toy when you put it in his

hand, uses her arm to swing at toys, brings hands to mouth, pushes up onto elbows/forearms

when on tummy

6 months Rolls from tummy to back, pushes up with straight arms when on tummy, leans on hands to

support himself when sitting

9 months Gets to a sitting position by herself, moves things from one hand to her other hand

¹Centers for disease control and prevention

12 Months



15 Months



18 Months



2 Years



12 months Pulls up to stand, walks, holding on to furniture

15 months Takes a few steps on his own


18 months Walks without holding on to anyone or anything, climbs on and off a couch or chair without help,

Kicks a ball, runs, walks (not climbs) up a few stairs with or without help 2 years

2) Mobility is a human right²

Delayed access to a mobility device enabling self-directed mobility restricts a child's right to mobility. Advocacy within the medical setting, with families & caregivers, and with funding sources is part of the process to obtain an appropriate mobility device for each child.

3) Four types of children can benefit from powered mobility³



Diagnostic examples


Spinal muscular atrophy types 1&2, congenital muscular dystrophy; multiple limb deficiencies; high level SCIs

Cerebral palsy, osteogenesis imperfecta, juvenile arthritis, myelomengeocele

Spina bifida, Duchennes muscular dystrophy, spinal muscular atrophy type 3, progressive neuromuscular


Cerebral palsy, Down Syndrome, orthopaedic issues earlier in life, ambulation delays

4) Propulsion Patterns & Best Practices⁴



Propulsion patterns

- A. SC-Semicircular
- B. SLOP-Single Loop Over C. DLOP- Double Loop Over
- D. ARC- Arc
- E. Resembles SLOP but hands raise above pushrim during recovery period.



²UN General Assembly, ³ Livingstone, R., & Paleg, G., ⁴ Slavens BA, Schnorenberg AJ, Aurit CM, Tarima S,

Propulsion pointers

- Pediatric propulsion differs from adults
- Children switch between stroke patterns during mobility and activities
- Consider all functional tasks when planning treatment and longer- term mobility strategies
- Wheelchair skill training has been shown to reduce fatigue and pain in children ages 8-18 using a modified Wheelchair skills Test⁵

5) Manual wheelchair configuration tips^{4,5,6}


- 1) Developmental differences between pediatrics & adults
 - a. Anatomy/neurological/growth
 - b. Physical endurance and fatigue
- 2) Seat child IN the wheelchair, not ON it
 - c. Age-appropriate seat-to-floor height
 - d. Center of Gravity (COG)
- 3) Ensure proper access to the wheels
 - e. Frame width (positioning, fit)
 - f. Rear seat-to-floor height (rear wheel access)
 - g. Camber (stability, rear wheel access)

- 4) Set realistic expectations
 - a. Age, development
 - b. Diagnosis
- 5) Avoid overuse of accessories
 - a. Consider weight of accessories
 - b. Consider necessity of accessories

- 6) Driving to Learn⁷: 8 Phases in 3 Stages (below)
 - 1. Explore **FUNCTIONS** body & tool/s
 - 2. Explore **SEQUENCING** body, tool/s & environment
 - 3. Explore **PERFORMANCE** body tool/s, environment, & occupation

7) Pediatric Power Wheelchair Training Best Practices^{8,9}

- Keep an open mind
- Age & IQ are not always good indicators to operate a power chair
- Avoid only focusing on readiness skills and pass-fail tests.



Occurs over time Should involve all key players Highly emotional process Patience, perseverance, & determination

Do¹⁰

- Encourage child to come to you/ towards something
- Provide positive & safe feedback
- Allow child time to figure out a solution rather than immediately providing solutions

Don't10

- · Instruct using directional commands
- Provide negative feedback
- Expect immediate mastery using powered mobility device

⁴Slavens BA, Schnorenberg AJ, Aurit CM, Tarima S, Vogel LC, Harris GF, ⁵Rammer, Jacob R, ⁶Boninger M., Souza A., Cooper R., Fitzgerald S., Koontz A., Fay B., ⁷Nilsson, L., & Durkin, J., ⁸ Bray, N., Kolehmainen, N., McAnuff, J., Tanner, L., Tuersley, L., Beyer, F., Craig, D., ⁹Kenyon, L. K., Mortenson, W. B., & Miller, W. C., ¹⁰Jones, M. A., McEwen, I. R., & Neas, B.R.

8) Justification

Letters of Medical Necessity (LMNs)

Proactively address:

- a. Child's age
- b. Child's capability of independence
- c. Why dependent mobility is not appropriate
- d. Age-appropriate supervision (child is a minor and need 24 hour supervision while engaging in mobility as do their peers).

Appealing LMN denials:

- a. Request denial in writing & clarify specific reason(s) for denial
- b. Address only denied items
- c. Include research to support your justification
- d. Be mindful of deadlines
- e. Fair hearings in-person and/or

Key Take Aways

- 1. Mobility is a human right
- 2. Benefits of "on-time" self-directed mobility
- 3. Strategies for enhancing mobility success in pediatrics
- 4. Funding & advocacy are part of the process

Suggested Next Steps

CEUs & workshops

Clinical Implications for Wheelchair Provision: Manual and Power Mobility

Wheelchair Configuration: The Importance of an Optimized Ride

TiLite Pilot, Twist, Aero X & SmartDrive

Explorer Mini

Permobil power wheelchairs- Koala, K450 MX, M300 PSJr, K300 PSJr

Comfort Inception cushion

Permobil Power Standing Training



Resources

Centers for Disease Control Developmental Milestones 2022

A Parent's Guide to the Explorer Mini

Clinician's Guide to Explorer Mini

The Wheelchair Handbook: A consumer's guide to seating & mobility equipment

A Guideline for Introducing Powered Mobility to Infants and Toddlers

RESNA Position on the Application of Power Mobility Devices for Pediatric Users- Update 2017

Alternate Funding Resources

Permobil LMN Generator

Permobil Clinical Resources and Clinical Research

References

- Boninger M., Souza A., Cooper R., Fitzgerald S., Koontz A., Fay B. (2002). Propulsion patterns and pushrim biomechanics in manual wheelchair propulsion. *Arch. Phys. Med. Rehabil.* 83, 718–723. 10.1053/apmr.2002.32455
- Bray, N., Kolehmainen, N., McAnuff, J., Tanner, L., Tuersley, L., Beyer, F., Craig, D. (2020). Powered mobility interventions for very
 young children with mobility limitations to aid participation and positive development: the EMPoWER evidence synthesis. *Health*Technology Assessment (Winchester, England), 24(50), 1.
- Centers for disease control and prevention. (n.d.). CDC's developmental milestone checklist. https://www.cdc.gov/ncbddd/actearly/milestones/index.html
- Jones, M. A., McEwen, I. R., & Neas, B.R. (2012). Effects of power wheelchairs on the development and function of young children with severe motor impairments. *Pediatr Phys Ther, 24*(2),131-140. doi: 10.1097/PEP.0b013e31824c5fdc.
- Kenyon, L. K., Mortenson, W. B., & Miller, W. C. (2018). 'Power in mobility': parent and therapist perspectives of the experiences of children learning to use powered mobility. *Developmental Medicine & Child Neurology, 60*(10), 1012-1017
- Livingstone, R., & Paleg, G. (2014). Practice considerations for the introduction and use of power mobility for children. *Developmental Medicine & Child Neurology*, 56(3), 210-221.
- Nilsson, L., & Durkin, J. (2014). Assessment of learning powered mobility use—Applying grounded theory to occupational performance. Journal of Rehabilitation Research and Development, 51(6), 963-974.
- Rammer, Jacob R. PhD; Krzak, Joseph J. PT, PCS, PhD; Slavens, Brooke A. PhD; Winters, Jack M. PhD; Riedel, Susan A. SM, PE; Harris, Gerald F. PhD, PE. Considering Propulsion Pattern in Therapeutic Outcomes for Children Who Use Manual Wheelchairs. Pediatric Physical Therapy: October 2019 Volume 31 Issue 4 p 360-368 doi: 10.1097/PEP.00000000000000649
- Slavens BA, Schnorenberg AJ, Aurit CM, Tarima S, Vogel LC, Harris GF. Biomechanics of Pediatric Manual Wheelchair Mobility. Front Bioeng Biotechnol. 2015 Sep 10;3:137. doi: 10.3389/fbioe.2015.00137. PMID: 26442251; PMCID: PMC4564732.
- UN General Assembly, Convention on the Rights of Persons with Disabilities: resolution / adopted by the General Assembly, 24 January 2007, A/RES/61/106, available at: https://www.refworld.org/docid/45f973632.html [accessed 20 December 2022]