






MAXIMIZING PROPULSION EFFICIENCY IN MANUAL WHEELCHAIRS: PRACTICAL STRATEGIES

CEU SIGN IN

PRESENTED BY: ERIN MANIACI, PT, DPT, ATP
ALPINE | OCTOBER 4, 2024




1

FACULTY DISCLOSURE

ERIN MANIACI, PT, DPT, ATP

- Physical Therapist
 - Inpatient SCI rehab, Amputee inpatient rehab, wheelchair prescription writing, gait training, outpatient neuro
- Current Clinical Education Specialist, Motion Composites
- Based in Phoenix, Arizona
- Grew up in MO → MU Grad → GO CHIEFS & TIGERS
- e.maniaci@motioncomposites.com



CEU SIGN IN

2




GREETINGS FROM MOTION COMPOSITES INDUSTRY AWARDS

CEU SIGN IN











3




CEU SIGN IN

MOTION COMPOSITES MISSION, VISION & VALUES

“LEAD THE EVOLUTION OF MOBILITY FOR BETTER LIVING”

“BECOME EVERYONE’S FAVORITE MOBILITY COMPANY BY OFFERING INNOVATIVE PRODUCTS AND THE MOST RESPONSIVE SERVICE”



4

Motion U

EVIDENCE BASED PRACTICE

OUR FOUNDATION

CLINICAL EXPERTISE

BEST RESEARCH EVIDENCE

PATIENT VALUES

EBP

Integrating individual clinical expertise with the best available external clinical evidence from systematic research.
-David Sacket, 1996

CEU SIGN IN

5

Motion U

LEARNING OBJECTIVES

BY THE END OF THIS COURSE PARTICIPANTS WILL BE ABLE TO:

1. Differentiate three aspects of manual wheelchair design and configuration that will influence propulsion efficiency.
2. Analyze three key points of manual wheelchair configuration their impact propulsion efficiency.
3. Implement two techniques to assess manual wheelchair propulsion on level surfaces.

6

Motion U

WHAT IS PROPULSION EFFICIENCY?

- Propulsion
 - the action of driving or pushing forward
- Efficiency
 - (especially of a system or machine) achieving maximum productivity with minimum wasted effort or expense

7

Motion U

Understanding Upper Extremity Injury

What the evidence tells us:

- Shoulder pain prevalence in manual wheelchair users can be upwards of 70-100%
- Nearly 75% of paraplegics with shoulder pain were suffering from impingement syndrome
- Degree of pain/derangement/disability is directly proportional to age/time since injury
- Shoulders are not designed to be weight-bearing joints

Guidelines/Recommendations

8

Motion U

GUIDELINES & RECOMMENDATIONS

- The Consortium of Spinal Cord Medicine (2005): "Preservation of Upper Limb Function Following Spinal Cord Injury: A Clinical Practice Guideline for Healthcare Professionals."
- RESNA Position Paper (2022): The Application of Ultralight Manual Wheelchairs
- Sawatzky, B. et. al. "The Need for Updated Clinical Practice Guidelines for Preservation of Upper Extremities in Manual Wheelchair Users: A Position Paper." American Journal of Physical Medicine and Rehabilitation, 2014.
- Cowan, Rachel E. Et.al. Impact of surface type, wheelchair weight and Axle Position on Wheelchair Propulsion by Novice Older Adults. Archives, Physical Medicine and Rehabilitation, 2009 July. 90(7): 1076-1083.

9

Motion U

THE CLINICAL GOAL

- Maximize Function
- Maximize/Improve Quality of Life (QOL)
- Minimize fatigue during normal daily activities
- Delayed onset or prevention of UE injury/dysfunction



10




CONSIDERATIONS FOR IMPROVING EFFICIENCY

MANUAL WHEELCHAIR SELECTION & TRAINING

11


Motion U

START WITH FRAME BASICS



FOLDING

- What most think of when they hear wheelchair, most common
- Typically, more weight
- More moving parts, less efficient
- Statistically chair of choice




RIGID

- Typically, lighter weight
- Rigidity improves efficiency
- More built-in adjustability
- More responsive to user input
- Easy to transport*

12

FOLDING FRAME STYLE HAS AN IMPACT


- Rigidity = less flex/movement in a frame which is more energy efficient
- Components of the folding frame design:
 - 2 x side frames
 - 1 x cross brace
 - Attachment of the cross brace
 - Seat saddles
- What details do we want to look for?



13

SIDE FRAMES

- Reduce potential for movement
- Less moving parts/bolts that loosen over time
- Lighter weight



14

CROSS BRACE


- Distributes weight evenly over the seat
- Reduces or eliminates torsion in the frame (more rigid)



15

OTHER FRAME COMPONENTS

- Oversized cross brace axles
 - More surface area
 - Stiffens ride; increases frame stability
- Locking seat saddles
 - Increases rigidity
 - Increases overall seat stability



16

Motion U

CHAIR MATERIALS

Typical materials: aluminum, titanium, carbon fiber

Material choices impact:

- Weight
 - Impacts forces at UE during startup/stop & turning
 - May limit functional transport
- Frame flexion/stiffness
 - Increased frame movement means energy is lost
 - Ride characteristics are impacted by materials ability to impact vibration transmitted to user

Impact of material is based on material science and the construction of the frame



17

Motion U

CONSIDER OVERALL WEIGHT

Standard definitions have not kept up with technology

US Medicare Definition of Ultralightweight = chair that weighs under 30 lbs. without wheels and has an adjustable COG

Folding vs. Rigid: impact of the cross brace

Every option weighs something:

- Armrests and brackets
- Anti-tips and brackets
- Push handles

18



PARAMETERS OF WHEELCHAIR CONFIGURATION

ASSUMING APPROPRIATE FRAME DIMENSIONS

19

Motion U

IMPORTANCE OF REAR WHEEL POSITION

THE EVIDENCE

Research shows that vertical and horizontal wheel position one of the most important adjustments to minimize impact on the UE during propulsion

- Medola, FO. Et. Al. Aspects of Manual Wheelchair Configuration Affecting Mobility: A Review. *Journal of Physical Therapy Science*. 26: 313-318, 2014.
- Rehabilitation Engineering & Assistive Technology Society of North America (RESNA). 2012. Position on the Application of Ultralight Manual Wheelchairs [position paper]. Retrieved from: RESNA: www.resna.org/resources/position_papers.dot
- Freixes O, Fernandez SA Gatti MA., Crespo, MJ, Olmos LE & Rubel IF (2010). Wheelchair axle position effect on start-up propulsion performance of persons with tetraplegia. *Journal of Rehabilitation Research & Development*, 47(7): 661-668.

20

Motion U

HORIZONTAL AXLE POSITION

- Move the rear axle as far forward as possible without compromising stability of the user
 - At or in front of the shoulder
 - Consider use of anti tips
 - Train clients in wheelies & fall recovery
- Axle position is USER Specific

REARWARD HORIZONTAL AXLE POSITION OPTIMAL HORIZONTAL AXLE POSITION

21

Motion U

FORWARD AXLE POSITION

22

Motion U

VERTICAL POSITION CHANGES

EXCESSIVELY LOW AXLE POSITION EXCESSIVELY HIGH AXLE POSITION

■ Path of Propulsion ■ Path of Propulsion ■ Shoulder Trauma

23

Motion U

VERTICAL AXLE POSITION

- Angle between the arm and forearm between 100-120 degrees (60-80 degrees of elbow flexion) with hand on top of pushrim
- Middle finger touches the center of the axle with arms hanging
- Can be difficult to optimize with some clients
- Look for:
 - Limited shoulder elevation
 - Reduced extension
 - Decreased external rotation


POSITION TEST 2 POSITION TEST 1

■ Path of Propulsion

24

Motion U

OPTIMAL REAR POSITION




- Results in
 - More efficient propulsion
 - The chair appears to roll easier and “feels lighter”

25

Motion U

REAR WHEEL SELECTION



Weight matters!
Stiffness = increased lateral stability; more efficiency?
Look at wheel movement during turns

Maintenance realities
Spoke wheels require little to no maintenance (user dependent)
Mag wheels warp over time

26

Motion U

REAR TIRE SELECTION




- The evidence is clear: air filled tires are more efficient
 - At 50% of recommended inflation, an air-filled tire rolls with less resistance than a solid tire. (Sawatzky et. al.)
- The maintenance reality
 - Tires should be checked weekly
 - Inflation needed 1-2x per month
 - New tire technology is puncture resistant
- Educate your client
 - Give them the responsibility & the choice**

27

Motion U

OTHER CONSIDERATIONS




- Front wheel (caster) selection
- Camber
- Handrims
- Seating and positioning options


28



29



30



31



32

Motion U

IDENTIFY PROPULSION PATTERNS

- Identify the current pattern
- Give instructions on most efficient movements
- Practice/re- evaluate
- QUESTION:** Which propulsion method do you think is the most efficient for traveling across flat terrain?

Semicircular Arc

Single Loop Over Double Loop Over

33

Motion U

SINGLE LOOPING OVER PROPULSION

Single Loop Over

Single looping over propulsion: the hands rise above the hand rim during recovery phase

34

Motion U

DOUBLE LOOP OVER

Double Loop Over

Double looping over propulsion: the hands rise above the hand rim, then cross over and drop below the hand rim during the recovery phase.

35

Motion U

ARCING


Arc

Arcing: The third metacarpophalangeal (MP) follows an arc along the path of the hand rim during the recovery phase

36

Motion U

SEMICIRCULAR



Semicircular: the hands fall below the hand rim during recovery phase


Semicircular

37

Motion U

PROTOCOL

- 10 m on level surfaces
- 10 m on carpet
 - Count number of strokes
 - Determine velocity
 - Observe pattern
- Up ADA ramp
- Figure 8
- Wheelies



38

Motion U

WHAT YOU NEED TO KNOW

You do have to teach the most efficient pattern
 Don't assume they know the best way
 Practice yourself
 Be creative; people learn differently
 Experienced users will self select the most advantageous pattern in difficult situations (i.e. ramps)
 Consider a propulsion training in plan of care

39

Motion U

WHEELCHAIR PROPULSION TEST

Wheelchair Propulsion Test (WPT) Version 1.0 Form

Subject #: _____ Date: _____ Time: _____ Test #: _____

Recorded Data*	
1. Able to successfully complete the 10m distance?	Yes <input type="checkbox"/> No <input type="checkbox"/>
2. Direction of travel	Forward <input type="checkbox"/> Backward <input type="checkbox"/>
3. Limbs contributing to propulsion, steering or braking (tick all that apply)	Left: Hand <input type="checkbox"/> Leg <input type="checkbox"/> Right: Hand <input type="checkbox"/> Leg <input type="checkbox"/>
4. Limb monitored for timing propulsion cycles (tick one limb)	Left: Hand <input type="checkbox"/> Leg <input type="checkbox"/> Right: Hand <input type="checkbox"/> Leg <input type="checkbox"/>
5. Time (to nearest second)	_____ s
6. Total number of propulsive cycles (to nearest full cycle)	_____ cycles
7. If using one or more hands for propulsion in the forward direction, during the correct phase, did the subject generally begin the contact between the hands and the hand-rim behind the up-front center of the test wheel?	Yes <input type="checkbox"/> No <input type="checkbox"/> Not applicable <input type="checkbox"/>
8. If using one or more hands for propulsion in the forward direction, during the recovery phase, did the subject generally use a path of the hands that was predominantly beneath the hand-rim?	Yes <input type="checkbox"/> No <input type="checkbox"/> Not applicable <input type="checkbox"/>
9. If using one or more feet for propulsion and going forward, did the subject make initial foot contact with the knee flexed less than 90° from full extension and finish with the knee flexed more than 90° (or the opposite if going backward)?	Yes <input type="checkbox"/> No <input type="checkbox"/> Not applicable <input type="checkbox"/>
10. Comments: (e.g., position on seat, trunk and arm posture, hand grip, foot contact, consistency, need for training, footwear, equipment worn, wheelchair issues)	_____

Derived Wheelchair-Propulsion Data*	
1. Speed: 10m / _____ # seconds =	_____ m/s
2. Push frequency (cadence): _____ # cycles / _____ # seconds =	_____ cycles/s
3. Effectiveness: 10m / _____ # cycles =	_____ m/cycle

*Directions on next page.

Tester signature: _____ Tester name (print): _____

40



**INTEGRATING INTO
CLINICAL PRACTICE**
BACK TO REALITY

41

Motion U

**CONSIDER CLIENT
SPECIFIC INDICATIONS**

- History of UE pain/dysfunction
- Reports of fatigue during the day
- Decline in independent mobility
- Diagnosis driven changes in function
- Caregiver abilities and needs*

42

Motion U

MAKING SENSE OF IT ALL

- Consider frame style carefully based on functional goals and abilities
- Pay attention to configuration details
- Integrate propulsion assessment/training into wheelchair delivery process
- Identify indications for follow up or further clinical intervention
- Provide your clients with the best opportunity for positive outcomes!


43

Motion U

**KEYS TO OPTIMIZING
PROPULSION**

1. **FRAME SELECTION**
2. **PROPULSION TRAINING**
3. **CONFIGURATION**


44



Thank You!

For more information contact:
Education@MotionComposites.com


46



RESOURCES

- Beekman CE, Miller-Porter L, Schoneberger M. Energy cost of propulsion in standard and ultralight wheelchairs in people with spinal cord injuries. *Phys Ther.* 1999;79(2):146-58.
- Boninger ML, Koontz AM, Sisto SA, Dyson-Hudson TA, Chang M, Price R et al. Pushrim biomechanics and injury prevention in spinal cord injury: Recommendations based on CULP-SCI investigations. *J Rehabil Res Dev* 2005;42(3 Suppl 1):9-20.
- Bossuyt, F. M., Boninger, M. L., Cools, A., Hogaboom, N., Eriks-Hoogland, I., & Arnet, U. (2020). Changes in supraspinatus and biceps tendon thickness: influence of fatiguing propulsion in wheelchair users with spinal cord injury. *Spinal Cord*, 58(3), 324-333. <https://doi.org/10.1038/s41393-019-0376-z>
- Caspall, J. J., Seligsohn, E., Dao, P. V., & Sprigle, S. (2013). Changes in inertia and effect on turning effort across different wheelchair configurations. *Journal of Rehabilitation Research and Development*, 50(10), 1353-1362. doi:10.1682/jrrd.2012.12.0219
- Cowan RE, Nash MS, Collinger JL, Koontz AM & Boninger ML (2009). Impact of surface type, wheelchair weight and axle position on wheelchair propulsion by novice older adults. *Archives of Physical Medicine and Rehabilitation*, 90(7): 1076-1083.
- Freixes O, Fernandez SA Gatti MA., Crespo, MJ, Olmos LE & Rubel IF (2010). Wheelchair axle position effect on start-up propulsion performance of persons with tetraplegia. *Journal of Rehabilitation Research & Development*, 47(7): 661-668.
- Louis, N., & Gorce, P. (2010). Surface electromyography activity of upper limb muscle during wheelchair propulsion: Influence of wheelchair configuration. *Clinical Biomechanics*, 25(9), 879-885. doi:10.1016/j.clinbiomech.2010.07.002

47



RESOURCES

- MacGillivray, M. K., Eng, J. J., Dean, E., & Sawatzky, B. J. (2020). Effects of Motor Skill-Based Training on Wheelchair Propulsion Biomechanics in Older Adults: A Randomized Controlled Trial. *Archives of Physical Medicine and Rehabilitation*, 101(1), 1-10. <https://doi.org/10.1016/j.apmr.2019.07.017>
- Medola, F. O., Elui, V. M., Santana, C. D., & Fortulan, C. A. (2014). Aspects of Manual Wheelchair Configuration Affecting Mobility: A Review. *Journal of Physical Therapy Science*, 26(2), 313-318. doi:10.1589/jpts.26.313
- Paralysed Veterans of America Consortium for Spinal Cord Medicine. (2005). Preservation of upper limb function following spinal cord injury: A clinical guideline for health-care professionals. *Journal of Spinal Cord Medicine*, 28(5):434-470.
- Requejo P, Mulroy, Haubert, et al. Evidence-Based Strategies to Preserve Shoulder Function in Manual Wheelchair Users with Spinal Cord Injury. *Top Spinal Cord Injury Rehabilitation*. 2008; 13(4): 86-119.
- Slowik, J. S., Requejo, P. S., Mulroy, S. J., & Neptune, R. R. (2016). The influence of wheelchair propulsion hand pattern on upper extremity muscle power and stress. *Journal of Biomechanics*, 49(9), 1554-1561. <https://doi.org/10.1016/j.jbiomech.2016.03.033>
- Stanfill, C. J., & Jensen, J. L. (2017). Effect of wheelchair design on wheeled mobility and propulsion efficiency in less-resourced settings. *African Journal of Disability (Online)*, 6, 1-9. <https://doi.org/10.4102/ajod.v6i0.342>
- Symonds, A., Barbareschi, G., Taylor, S., & Holloway, C. (2018). A systematic review: the influence of real time feedback on wheelchair propulsion biomechanics. *Disability & Rehabilitation: Assistive Technology*, 13(1), 47-53. <https://doi.org/10.1080/17483107.2016.1278477>
- Sawatzky, Bonita & Denton, Ian. (2006). Wheeling efficiency: The effects of varying tyre pressure with children and adolescents. *Pediatric Rehabilitation*, 9, 122-6. doi:10.1080/13638490500126707
- Tafner Curi, H., de Lima, J., & Chaves Ferretti, E. (2020). Factors related to propulsion efficiency in manual wheelchair users with paraplegia due to spinal cord injury. *Brazilian Journal of Occupational Therapy / cadernos brasileiros de Terapia Ocupacional*, 28(3), 999-1019. <https://doi.org/10.4322/2576-8910.ctoAB1935>

48