



MANUAL WHEELCHAIR PROPULSION: TECHNIQUES, POSITIONING, AND EQUIPMENT CONSIDERATIONS

PRESENTED BY: Petra Conaway, PT, PT, ATP/SMS
DATE: OCTOBER 2026



FACULTY DISCLOSURE

PETRA CONAWAY, PT, DPT, ATP/SMS

- Clinical Education Specialist, Motion Composites
- Physical Therapist
 - Wheelchair & Seating Specialist @ CIDE, Denver, CO
 - Inpatient/Outpatient SCI rehab @ Craig Hospital, Denver, CO
 - Research in SCI, CVA and CP @ RIC, Chicago, IL
- Czech Republic → England → Illinois → Colorado

p.conaway@motioncomposites.com





MOTION COMPOSITES



“LEAD THE EVOLUTION OF
MOBILITY FOR BETTER LIVING”




 *Wishes*
for Wheels



EVIDENCE BASED PRACTICE

OUR FOUNDATION



Integrating individual clinical expertise with the best available external clinical evidence from systematic research.

-David Sacket, 1996

LEARNING OBJECTIVES

BY THE END OF THIS COURSE, PARTICIPANTS WILL BE ABLE TO:

1. Analyze the impact of client ability, musculoskeletal health, and long-term mobility outcomes on the efficiency of different propulsion techniques.
2. Evaluate client positioning and wheelchair configuration to enhance propulsion mechanics, minimize potential for injury, and maximize independence.
3. Apply principles of component selection to customize wheelchair setup to support injury prevention and long-term functional mobility with propulsion.





WHAT IS PROPULSION

?





WHAT IS PROPULSION

- The method a person uses to move their wheelchair forward, backward, and during turns
- Can be done using arms/hands, feet/legs, or a combination
- Efficiency and technique affect speed, energy use, comfort, and long-term health (e.g., shoulders, posture)
- GOAL: match the wheelchair and propulsion method to the individual's abilities, environment, and daily activities





TYPES OF PROPULSION

- Arm
- Foot
- Hemi
- Dependent



A grayscale close-up of a bicycle wheel, showing the spokes and the tire. A red horizontal bar is positioned above the text.

ARM PROPULSION



ARM PROPULSION

- Using arms and hands typically by pushing on the handrim
 - Most common method of manual wheelchair propulsion
- Energy is applied to the system when arms are in contact with the handrim/tire
 - **How can you be efficient?**
 - Increase the contact time – hand on handrim
- Clients can push on their handrim 2000-3000 times a day!
 - Not just young & active ones
- Consideration: High risk of shoulder pain or overuse injury if not optimized
 - Shoulders are not meant to be weight bearing joints



RESEARCH & LITERATURE


- High demand on shoulders, wrists, and hands
- Upper extremity pain and strain affects up to 76% of manual wheelchair users
- Risk of overuse injury
 - Shoulder, wrist, elbows
 - Rotator cuff degeneration, shoulder impingement syndrome, tendonitis and bursitis, carpal tunnel syndrome
- Most research is on people with SCI






RESEARCH & LITERATURE

- A lighter wheelchair *can* improve efficiency and reduce strain, but the optimal configuration includes a combination of ideals across multiple aspects.
 - Lightest weight
 - Optimal set-up
 - High-quality components
 - Propulsion technique





PHYSICAL CONSIDERATIONS FOR ARM PROPULSION





OPTIMAL POSITIONING


- Rear wheel set up
 - Dictates arm & elbow positioning
 - Allows better contact with more of the handrim
 - Important for all populations
 - Don't forget the geriatric client
- Two ways to achieve this optimal positioning
 - Rear wheel size
 - Positioning of wheels within the frame






HORIZONTAL REAR AXLE POSITION


- Horizontal axle position = Center of gravity (COG)
- Move the rear wheel as far forward as possible without compromising the stability of the user
- Axle at or in front of the shoulder
- Axle position is USER Specific




REARWARD HORIZONTAL AXLE POSITION



OPTIMAL HORIZONTAL AXLE POSITION






FORWARD AXLE POSITION

- Allows access to the full wheels for long active phase on push rim
- Allows for use of multiple muscle groups for endurance and power




VERTICAL REAR AXLE POSITION


- Angle between the arm and forearm between 100-120°
 - 60-80° of elbow flexion when the hand is on the top of the pushrim
- Middle finger touches the center of the axle with arms hanging
- Can be difficult to optimize with some clients



EXCESSIVELY HIGH AXLE POSITION



EXCESSIVELY LOW AXLE POSITION



IDEAL POSITION




PRACTICAL APPLICATION


- LIMITATIONS TO ACHIEVING OPTIMAL POSITIONING
 - Transfer height
 - Environment
 - Environmental surface heights too low (tables, counters, etc.)
 - Clearance (knee & ground)
 - Seating surface heights
 - Safety with forward COG
- TROUBLE SHOOTING
 - Modify seat to floor height
 - Rear or both front & rear
 - Modify rear wheel size

IMPACT OF REAR SEAT TO FLOOR HEIGHT

1.5" DIFFERENCE



3" DIFFERENCE



IMPACT OF OPTIMAL CONFIGURATION





WHEELCHAIR MODEL CONSIDERATIONS

FRAME HEIGHT

- Can you achieve optimal STFH
- Consider future adjustment up or down

VERTICAL ADJUSTABILITY

- Can you achieve optimal STFH
- Increments for modification: 1", 0.5", ¼", infinite

HORIZONTAL ADJUSTABILITY

- As forward as possible without compromising safety
- How far forward can it go
- Increments for modification: 1", 0.5", ¼", infinite

CASTERS & REAR WHEELS

- Will dictate STF heights
- Will impact maneuverability
- Rear wheel: impact arm positioning

SEATING SYSTEM

- Seat cushion: Height will impact positioning in the frame and user interacting in the environment
- Back support: should not impact shoulder ROM, PSIS support

HANDRIM SELECTION

- Material
- Shape
- Size
- Spacing



Semicircular



Arc



Single Loop
Over



Double Loop
Over

TECHNIQUE: PROPULSION PATTERNS

- 4 main upper extremity propulsion patterns
- Vary in efficiency depending on the situation
 - Functional Goal
 - Environment



SINGLE LOOP OVER



ACTION

Hands grab the handrim, push down and then rise above the hand rim during recovery phase

PRACTICAL APPLICATION

Crowded environments

RISK

Shoulder over use with elevation

SINGLE LOOP OVER PROPULSION





DOUBLE LOOP OVER



ACTION

Hands grab the handrim, push down, then let go with wrist flicking and drop below the hand rim during the recovery phase

PRACTICAL APPLICATION

Fast starts, when on rough terrain, carpet etc.


RISK

Wrist overuse – carpal tunnel

DOUBLE LOOP OVER PROPULSION



ARCING



ACTION

The hand does not leave the handrim during this style leading to a slow down during the recovery phase

PRACTICAL APPLICATION

Inclines

RISK

Overuse injury due to inefficiency, most inherent pattern people will choose – risk of never getting more efficient

ARCING PROPULSION



SEMICIRCULAR



ACTION

Hands grab the back of the handrim, push forward and drop below the hand rim during recovery phase

PRACTICAL APPLICATION

Most efficient for everyday use

RISK

Minimized on all aspects

SEMICIRCULAR PROPULSION



A grayscale photograph of a bicycle wheel, showing the spokes and the tire. A red horizontal bar is positioned above the text.

FOOT PROPULSION



FOOT PROPULSION

- Moving a wheelchair with the lower extremities using one or both feet
- Energy is applied to the system when foot is in contact with the ground
- Uses larger muscle groups
 - Typical weight bearing joints
 - Removes the high risk of upper extremity strain
 - More efficient and less effortful than hand propulsion
- Optimal postural alignment and positioning is just as important for foot propulsion

<https://www.youtube.com/watch?v=rKy7-npLXAg>



POPULATION

Common among older adults

Clients with limited UE strength or mobility

Previously ambulatory clients

- Prefer foot propulsion due to familiarity with using their LEs for movement

Clients with good hip/knee ROM, ankle mobility, and gluteal strength






RESEARCH & LITERATURE

- Under-researched topic
- Most studies use able-bodied subjects
 - Due to co-morbidities
- Lowering seat to floor height improves propulsion speed and reduces effort
- Anterior tilt (3-6°) facilitates more effective foot propulsion, hip & knee flexion ROM and user comfort






PHYSICAL CONSIDERATIONS OF FOOT PROPULSION





IMPACT OF THE LOWER EXTREMITIES

- Strongly influence pelvic mobility
- Example: Tight hamstrings
 - Affect the amount of isolated hip flexion
 - May promote sacral sitting
- MUST be considered for back support selection





OPTIMAL POSITIONING

- SEAT DEPTH
 - Maximize posterior thigh support
 - Allow for optimal power & recovery phases
 - Posterior buttocks to popliteal fossa minus 1-2"
- SEAT TO FLOOR HEIGHT (STFH)
 - Lowering STFH improves propulsion speed and reduces effort
 - Finished front STFH 1-2" lower than lower leg length
 - Compensates for the extra distance needed to achieve heel strike





OPTIMAL POSITIONING

- ANTERIOR SEAT SLOPE
 - 3-6°
 - ➔ 1-2" on an 18" seat depth
 - Pair with optimal support surface
 - Reduce risk of posterior pelvic tilt
 - Minimize increased pressure on coccyx & ITs

$$\text{Seat Slope (degrees)} = \arctan \left(\frac{Hr - Hf}{D} \right)$$



POSTERIOR V. ANTERIOR SEAT SLOPE

Posterior Seat Slope 	Anterior Seat Slope 
• Elevates knees and shifts pelvis rearward • Reduces foot contact on the ground	• Facilitates better foot/ground contact
• Increases risk of forward sliding	• Improves postural alignment
• Limits hip and knee ROM	• Enhances propulsion mechanics
• Compromises postural stability • Posterior pelvic tilt → kyphosis	• Reduces forward sliding and shear
• Reduces propulsion efficiency • Increased energy expenditure	



PRACTICAL APPLICATION



- LIMITATIONS TO ACHIEVING OPTIMAL POSITIONING
 - Transfer height
 - Environment
 - Environmental surface too high (tables, counters, etc.)
 - Ground clearance of footplates
 - Positioning with footplates
 - Seating surface heights



PRACTICAL APPLICATION

- TROUBLE SHOOTING
 - Center of Gravity
 - Common mistake is placing COG in a rearward position
 - Arms not being used
 - Age biased
 - Increases weight on the casters decreasing efficiency
 - Modify rear wheel size as needed
 - Larger diameter wheel facilitates maneuverability
 - Smaller diameter wheel can accommodate lower STFHS
 - Footwear
 - Backwards
 - Larger muscles used for power phase





WHEELCHAIR MODEL CONSIDERATIONS

FRAME

- Can you achieve optimal STFH
 - Consider future adjustment up or down
- Removeable hangers
- Will they be a dependent user in the future

VERTICAL ADJUSTABILITY

- Can you achieve optimal STFH
- Increments for modification
- 1", 0.5", ¼", infinite

HORIZONTAL ADJUSTABILITY

- How far forward can it go
- Increments for modification
- 1", 0.5", ¼", infinite

CASTERS & REAR WHEELS

- Will dictate STFHS
- Will impact maneuverability

SEATING SYSTEM

- Seat cushion: sufficient buttock well, pre-ischial shelf, cover considerations
- Back support: PSIS support



OTHER PROPULSION STYLES



HEMI PROPULSION



https://youtu.be/kQ3cvQNuCPQ?si=2qbBFDHI_1D3pEbt



<https://youtu.be/4WAlixzQRT4?si=ZVz5MKbVrsbS-Y4y>



HEMI PROPULSION



HEMI PROPULSION

WHAT IS IT?

- To propel using one hand and one foot, typically on one side of the body

POPULATION

- Often hemiplegia, paralysis or weakness on one side of the body

RESEARCH

- Similar movement patterns are seen in able-bodied and wheelchair users

OPTIMAL POSITIONING

- Utilizing principles from both arm and foot propulsion

LIMITATIONS

- There is always a compromise between upper and lower extremity set-up
- Propulsion efficiency is often slower and more fatiguing

TROUBLESHOOTING

- When compromising – favour the stronger extremity

WHEELCHAIR

- Removeable hangers
- One arm drive
- Wheel lock extension/Unilateral wheel-lock
- Pay attention to asymmetrical seating needs: Trunk/core fatigue is common due to high stabilization demand
- Handrim – because of only using one hand

TRAINING

- Apply both principles of arm and foot propulsion training



DEPENDENT PROPULSION



Central Grampians
Palliative Care

**Push
in a
wheelchair**



<https://www.youtube.com/watch?v=QleLHZ7djj0>



DEPENDENT PROPULSION

WHAT IS IT?

- When a caregiver/attendant propels the wheelchair

POPULATION

- User who is unable to achieve a successful or efficient propulsions stroke with upper or lower extremities

WHEELCHAIR


- Attendant wheel lock
- Wheel size can enhance or limit maneuverability
- Ensure optimal seated posture
- Center of Gravity
 - Efficient propulsion for attendant
 - Smaller wheelbase for maneuverability
- Appropriate height/style of push handles across all caregivers

TRAINING


- Assess caregiver(s) propelling user within the wheelchair

PROPULSION TRAINING

- **TEACH PROPULSION STYLE**
- Educate on risks of propulsion styles
- Assess overall WC Skills
 - Level surfaces
 - Up/down ramps
 - Carpet
 - Turning
 - Wheelies



WHAT YOU NEED TO KNOW



- You do have to teach the most efficient pattern
 - Just because they move does not mean it is the best way
 - Practice yourself
 - Be creative; people learn differently
- Experienced users will self-select the most advantageous pattern in difficult situations (i.e., ramps)
- Consider a propulsion training in plan of care



WHEELCHAIR PROPULSION TEST

Wheelchair Propulsion Test (WPT)[®] Version 1.0 Form

Subject #: _____ Date: _____ Time: _____ Test # _____

Recorded Data*	
1. Able to successfully complete the 10m distance?	Yes ☐ No ☐
2. Direction of travel	Forward ☐ Backward ☐
3. Limbs contributing to propulsion, steering or braking (tick all that apply)	Left: Hand ☐ Leg ☐ Right: Hand ☐ Leg ☐
4. Limb monitored for timing propulsion cycles (tick one limb)	Left: Hand ☐ Leg ☐ Right: Hand ☐ Leg ☐
5. Time (to nearest second)	_____ s
6. Total number of propulsive cycles (to nearest full cycle)	_____ cycles
7. If using one or more hands for propulsion in the forward direction, during the <i>contact phases</i> , did the subject generally begin the contact between the hands and the hand-rims behind the top dead center of the rear wheel?	Yes ☐ No ☐ Not applicable ☐
8. If using one or more hands for propulsion in the forward direction, during the <i>recovery phases</i> , did the subject generally use a path of the hands that was predominantly beneath the hand-rims?	Yes ☐ No ☐ Not applicable ☐
9. If using one or more <i>feet for propulsion</i> and going forward, did the subject make initial foot contact with the knee flexed less than 90° from full extension and finish with the knee flexed more than 90° (or the opposite if going backward)?	Yes ☐ No ☐ Not applicable ☐
10. Comments: (e.g., position on seat, trunk and arm posture, hand grip, foot contact, consistency, need for training, footwear, equipment worn, wheelchair issues)	
Derived Wheelchair-Propulsion Data*	
1. Speed: 10m / _____ # seconds =	_____ m/s
2. Push frequency (cadence): _____ # cycles / _____ # seconds =	_____ cycles/s
3. Effectiveness: 10m / _____ # cycles =	_____ m/cycle

*Directions on next page.

Tester signature: _____ Tester name (print): _____

Derived Wheelchair-Propulsion Data*

1. Speed: 10m / _____ # seconds =	_____ m/s
2. Push frequency (cadence): _____ # cycles / _____ # seconds =	_____ cycles/s
3. Effectiveness: 10m / _____ # cycles =	_____ m/cycle

*Directions on next page.

Tester signature: _____ Tester name (print): _____

*Propulsion Test*



MAKING SENSE OF IT ALL

- Consider frame style carefully based on functional goals and abilities
- Pay attention to configuration details
- Integrate propulsion assessment & training into wheelchair service provision
- Identify concerns for follow up or further clinical intervention
- Provide your clients with the best opportunity for positive outcomes!



Handouts





Thank You!

For more information contact:

Education@MotionComposites.com



REFERENCES

- Bossuyt, F. M., Boninger, M. L., Cools, A., Hogaboom, N., Eriks-Hoogland, I., & Arnet, U. (2020). Changes in supraspinatus and biceps tendon thickness: influence of fatiguing propulsion in wheelchair users with spinal cord injury. *Spinal Cord*, 58(3), 324–333. <https://doi.org/10.1038/s41393-019-0376-z>
- Charbonneau R, Kirby RL, Thompson K. (2013). Manual wheelchair propulsion by people with hemiplegia: within-participant comparisons of forward versus backward techniques. *Arch Phys Med Rehabil*. 2013 Sep;94(9):1707-13. doi: 10.1016/j.apmr.2013.03.001. Epub 2013 Mar 13. PMID: 23500180.
- Hamstra, C. , Jones, D. (2025). Considerations When Working with the Geriatric Population. In M. Lange & J. Minkel (Eds.) *Seating and Wheeled Mobility: A Clinical Resource Guide* (2nd ed., pp. 493-514). Routledge. <https://doi.org/10.4324/978100352637>
- Heinrichs, N. D., Kirby, R. L., Smith, C., Russell, K. F. J., Theriault, C. J., & Doucette, S. P. (2021). Effect of seat height on manual wheelchair foot propulsion, a repeated-measures crossover study: Part 1 — wheeling forward on a smooth level surface. *Disability and Rehabilitation: Assistive Technology*, 16(8), 831–839. <https://doi.org/10.1080/17483107.2020.1741036>
- Heinrichs, N. D., Kirby, R. L., Smith, C., Russell, K. F. J., Theriault, C. J., & Doucette, S. P. (2022). Effect of seat height on manual wheelchair foot propulsion, a repeated-measures crossover study: part 2 – wheeling backward on a soft surface. *Disability & Rehabilitation: Assistive Technology*, 17(3), 325–330. <https://doi.org/10.1080/17483107.2020.1782490>
- Lange, M. L., & Minkel, J. (2025) *Seating and Wheeled Mobility : A Clinical Resource Guide*. Routledge, 2025.
- MacGillivray, M. K., Eng, J. J., Dean, E., & Sawatzky, B. J. (2020). Effects of Motor Skill-Based Training on Wheelchair Propulsion Biomechanics in Older Adults: A Randomized Controlled Trial. *Archives of Physical Medicine and Rehabilitation*, 101(1), 1–10. <https://doi.org/10.1016/j.apmr.2019.07.017>

REFERENCES *continued*

- Tafner Curi, H., de Lima, J., & Chaves Ferretti, E. (2020). Factors related to propulsion efficiency in manual wheelchair users with paraplegia due to spinal cord injury. *Brazilian Journal of Occupational Therapy / Cadernos Brasileiros de Terapia Ocupacional*, 28(3), 999–1019. <https://doi.org/10.4322/2526-8910.ctoAR1935>
- da Silva Bertolaccini G, Sandnes FE, Medola FO, Gjøvaag T. (2022) Effect of Manual Wheelchair Type on Mobility Performance, Cardiorespiratory Responses, and Perceived Exertion. *Rehabil Res Pract*. 2022 Jun 11;2022:5554571. doi: 10.1155/2022/5554571. PMID: 35726217; PMCID: PMC9206577.
- Requejo PS, Furumasu J, Mulroy SJ. (2015). Evidence-Based Strategies for Preserving Mobility for Elderly and Aging Manual Wheelchair Users. *Top Geriatr Rehabil*. 2015 Jan-Mar;31(1):26-41. doi: 10.1097/TGR.0000000000000042. PMID: 26366040; PMCID: PMC4562294.
- Rehabilitation Engineering & Assistive Technology Society of North America [RESNA] (2022). Position on the Application of Ultralight Manual Wheelchairs [position paper]. Retrieved from: <https://www.resna.org/Portals/0/Position%20Papers/RESNA%20Position%20on%20the%20Application%20of%20Ultralight%20Manual%20Wheelchairs.pdf?ver=D1H8hccXxkWkG6D73Gh6WA%3D%3D>
- Soleymani, H., & Cowan, R. (2024). Relationship between rolling resistance, preferred speed, and manual wheelchair propulsion mechanics in non-disabled adults. *Disability and Rehabilitation: Assistive Technology*, 19(5), 1980–1991. <https://doi.org/10.1080/17483107.2023.2239297>